

StockOptions™ Salt is a kit preformulated, sterile filtered stock salt solutions for biological macromolecular crystallization. The StockOptions Salt reagents are supplied in varying concentrations, each concentration appropriate for each particular salt's application as a primary or secondary precipitant, or additive in a crystallization experiment. StockOptions Salt is comprised of 49 unique reagents.

### Suggested Use

StockOptions Salt is designed to help researchers improve the speed, accuracy, precision, and quality of the formulation of crystallization screen solutions and crystallization optimization solutions. Researchers can use the individual StockOptions reagents to conveniently formulate customer screen solutions or accurately reproduce standard screen solutions from Hampton Research kits such as Crystal Screen™, Crystal Screen Cryo™, Crystal Screen Lite™, Natrix™, Crystal Screen 2™, and PEG/Ion™. StockOptions Salt reagents can also be used to create solutions for the refinement and optimization of preliminary crystallization conditions. Finally, StockOptions Salt reagents can be used to create accurate, precise, reproducible, high quality solutions for the production of single crystals.

### Example 1

Crystal Screen Reagent 2 (1 ml plate reservoir)

Solution composition:

0.4 M Potassium sodium tartrate tetrahydrate.

Suggested Stock Solutions:

1.5 M Potassium sodium tartrate tetrahydrate.

- 1.) Pipet 267 microliters of 1.5 M Potassium sodium tartrate tetrahydrate into the plate reservoir.
- 2.) Pipet 733 microliters of sterile filtered deionized water into the plate reservoir.
- 3.) Aspirate and dispense the solution until homogenous.

### Example 2

Crystal Screen 2 Reagent 1 (1 ml plate reservoir)

Solution composition:

2.0 M Sodium chloride, 10% w/v PEG 6,000.

Suggested Stock Solutions:

5.0 M Sodium chloride, 50% w/v PEG 6,000.

- 1.) Pipet 400 microliters of sterile filtered deionized water into the plate reservoir.
- 2.) Pipet 200 microliters of 50% PEG 6,000 into the plate reservoir.
- 3.) Pipet 400 microliters of 5.0 M Sodium chloride into the plate reservoir.
- 4.) Aspirate and dispense the solutions until homogenous.

### Example 3

Crystal Screen Lite reagent 48 (1 ml plate reservoir)

Solution composition:

0.1 M TRIS hydrochloride pH 8.5, 1.0 M Ammonium phosphate monobasic

Suggested Stock Solutions:

StockOptions 1.0 M TRIS hydrochloride pH 8.5, 2.5 M Ammonium phosphate monobasic.

- 1.) Pipet 500 microliters of sterile filtered deionized water into the plate reservoir.
- 2.) Pipet 100 microliters of 1.0 M TRIS hydrochloride pH 8.5 into the plate reservoir.
- 3.) Pipet 400 microliters of 2.5 M Ammonium phosphate monobasic into the plate reservoir.
- 4.) Aspirate and dispense the solution until homogeneous.

### Example 4

A custom reagent of 1.7 M Sodium malonate, 0.1 M HEPES sodium pH 7.5.

Suggested Stock Solutions:

3.4 M Sodium malonate pH 7.0, StockOptions 1.0 M HEPES sodium pH 7.5.

- 1.) Pipet 400 microliters of sterile filtered deionized water into the plate reservoir.
- 2.) Pipet 100 microliters of 1.0 M HEPES sodium pH 7.5 into the plate reservoir.
- 3.) Pipet 500 microliters of 3.4 M Sodium malonate pH 7.0 into the plate reservoir.
- 4.) Aspirate and dispense the solution until homogeneous.

### For Best Results

Use Hampton Research StockOptions Salt together with StockOptions pH Screens and other Optimize™ reagents for the best results. The reagent in StockOptions Salt are available in larger volumes from the Hampton Research Optimize line of reagents should larger volumes be required.

StockOptions Salt reagents are stable at room temperature and are best if used within 12 months of receipt. To enhance reagent stability it is strongly recommended that StockOptions Salt be stored at 4°C or -20°C. Avoid ultraviolet light to preserve reagent stability.

### Technical Support

Inquiries regarding StockOptions Salt Screen reagent formulation, interpretation of screen results, optimization strategies and general inquiries regarding crystallization are welcome. Please e-mail, fax, or telephone your request to Hampton Research. Fax and e-mail Technical Support are available 24 hours a day. Telephone technical support is available 8:00 a.m. to 4:30 p.m. USA Pacific Standard Time.

Reagent Number	Reagent	Molecular Formula	Molecular Weight	CAS Number
1	1.0 M Ammonium acetate	CH <sub>3</sub> COONH <sub>4</sub>	77.08	[631-61-8]
2	5.0 M Ammonium chloride	NH <sub>4</sub> Cl	53.49	[12125-02-9]
3	2.5 M Ammonium phosphate monobasic	NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub>	115.03	[7722-76-1]
4	10.0 M Ammonium fluoride	NH <sub>4</sub> F	37.04	[12125-01-8]
5	10.0 M Ammonium formate	HCOONH <sub>4</sub>	63.06	[540-69-2]
6	2.5 M Ammonium citrate dibasic	(NH <sub>4</sub> ) <sub>2</sub> C <sub>6</sub> H <sub>6</sub> O <sub>7</sub>	226.19	[3012-65-5]
7	3.5 M Ammonium phosphate dibasic	(NH <sub>4</sub> ) <sub>2</sub> HPO <sub>4</sub>	132.06	[7783-28-0]
8	10.0 M Ammonium nitrate	NH <sub>4</sub> NO <sub>3</sub>	80.04	[6484-52-2]
9	3.5 M Ammonium sulfate	(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	132.14	[7783-20-2]
10	2.0 M Ammonium tartrate dibasic	(NH <sub>4</sub> ) <sub>2</sub> C <sub>4</sub> H <sub>4</sub> O <sub>6</sub>	184.15	[3164-29-2]
11	1.0 M Calcium acetate hydrate	(CH <sub>3</sub> COO) <sub>2</sub> Ca • xH <sub>2</sub> O	158.17	[114460-21-8]
12	2.0 M Calcium chloride dihydrate	Cl <sub>2</sub> Ca • 2H <sub>2</sub> O	147.01	[10035-04-8]
13	5.0 M Lithium acetate dihydrate	CH <sub>3</sub> COOLi • 2H <sub>2</sub> O	102.02	[6108-17-4]
14	10.0 M Lithium chloride	LiCl	42.39	[7447-41-8]
15	1.5 M Lithium citrate tribasic tetrahydrate	HOC(COOLi)(CH <sub>2</sub> COOLi) <sub>2</sub> • 4H <sub>2</sub> O	281.99	[6080-58-6]
16	8.0 M Lithium nitrate	LiNO <sub>3</sub>	68.95	[7790-69-4]
17	2.0 M Lithium sulfate monohydrate	Li <sub>2</sub> SO <sub>4</sub> • H <sub>2</sub> O	127.96	[10102-25-7]
18	1.0 M Magnesium acetate tetrahydrate	(CH <sub>3</sub> COO) <sub>2</sub> Mg • 4H <sub>2</sub> O	214.46	[16674-78-5]
19	2.0 M Magnesium chloride hexahydrate	MgCl <sub>2</sub> • 6H <sub>2</sub> O	203.30	[7791-18-6]
20	1.0 M Magnesium formate dihydrate	Mg(CHOO) <sub>2</sub> • 2H <sub>2</sub> O	150.38	[557-39-1]
21	3.0 M Magnesium nitrate hexahydrate	Mg(NO <sub>3</sub> ) <sub>2</sub> • 6H <sub>2</sub> O	256.41	[13446-18-9]
22	2.5 M Magnesium sulfate hydrate	MgSO <sub>4</sub> • aq	120.37	[22189-08-8]
23	4.0 M Nickel(II) chloride hexahydrate	NiCl <sub>2</sub> • 6H <sub>2</sub> O	237.69	[7791-20-0]
24	5.0 M Potassium acetate	CH <sub>3</sub> COOK	98.14	[127-08-2]
25	4.0 M Potassium chloride	KCl	74.55	[7447-40-7]
26	2.5 M Potassium citrate tribasic monohydrate	C <sub>6</sub> H <sub>5</sub> K <sub>3</sub> O <sub>7</sub> • H <sub>2</sub> O	324.41	[6100-05-6]
27	1.5 M Potassium phosphate monobasic	KH <sub>2</sub> PO <sub>4</sub>	136.09	[7778-77-0]
28	6.0 M Potassium fluoride	KF	58.10	[7789-23-3]
29	14.0 M Potassium formate	HCOOK	84.12	[590-29-4]
30	3.0 M Potassium phosphate dibasic	K <sub>2</sub> HPO <sub>4</sub>	174.18	[7758-11-4]
31	2.0 M Potassium nitrate	KNO <sub>3</sub>	101.10	[7757-79-1]
32	1.5 M Potassium sodium tartrate tetrahydrate	C <sub>4</sub> H <sub>4</sub> KNaO <sub>6</sub>	282.22	[6381-59-5]
33	0.5 M Potassium sulfate	K <sub>2</sub> SO <sub>4</sub>	174.26	[7778-80-5]
34	8.0 M Potassium thiocyanate	KSCN	97.18	[333-20-0]
35	3.0 M Sodium acetate trihydrate	CH <sub>3</sub> COONa • 3H <sub>2</sub> O	136.08	[6131-90-4]
36	5.0 M Sodium chloride	NaCl	58.44	[7647-14-5]
37	1.6 M Sodium citrate tribasic dihydrate	C <sub>6</sub> H <sub>5</sub> Na <sub>3</sub> O <sub>7</sub> • 2H <sub>2</sub> O	294.10	[6132-04-3]
38	5.0 M Sodium phosphate monobasic	NaH <sub>2</sub> PO <sub>4</sub> • H <sub>2</sub> O	137.99	[10049-21-5]
39	0.8 M Sodium fluoride	NaF	41.99	[7681-49-4]
40	7.0 M Sodium formate	HCOONa	68.01	[141-53-7]
41	1.0 M Sodium phosphate dibasic dihydrate	Na <sub>2</sub> HPO <sub>4</sub> • 2H <sub>2</sub> O	177.99	[10028-24-7]
42	3.4 M Sodium malonate pH 7.0	C <sub>3</sub> H <sub>4</sub> O <sub>4</sub> + NaOH	104.06 + NaOH	[141-82-2]
43	7.0 M Sodium nitrate	NaNO <sub>3</sub>	84.99	[7631-99-4]
44	1.0 M Sodium sulfate decahydrate	Na <sub>2</sub> SO <sub>4</sub> • 10H <sub>2</sub> O	322.20	[7727-73-3]
45	1.5 M Sodium tartrate dibasic dihydrate	C <sub>4</sub> H <sub>4</sub> Na <sub>2</sub> O <sub>6</sub> • 2H <sub>2</sub> O	230.08	[6106-24-7]
46	8.0 M Sodium thiocyanate	NaSCN	81.07	[540-72-7]
47	1.2 M Succinic acid pH 7.0	C <sub>4</sub> H <sub>6</sub> O <sub>4</sub> + NaOH	118.09 + NaOH	[110-15-6]
48	1.0 M Zinc acetate dihydrate	Zn(CH <sub>3</sub> COO) <sub>2</sub> • 2H <sub>2</sub> O	219.52	[5970-45-6]
49	2.0 M Zinc sulfate heptahydrate	ZnSO <sub>4</sub> • 7H <sub>2</sub> O	287.56	[7446-20-0]

Hampton Research  
34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

Technical Support e-mail: tech@hrmail.com

Website: www.hamptonresearch.com