

Teklad LM-485 Mouse/Rat Sterilizable Diet

Product Description- 7012 is a fixed formula, autoclavable diet manufactured with high quality ingredients and designed to support growth and reproduction of rodents. Typical isoflavone concentrations (daidzein + genistein aglycone equivalents) range from 300 to 600 mg/kg. Absence of animal protein and fish meal minimizes the presence of nitrosamines. 7012 is supplemented with additional vitamins to ensure nutritional adequacy after autoclaving. **Also available certified (7012C) and irradiated (7912).**

Ingredients (in descending order of inclusion)- Ground corn, dehulled soybean meal, ground oats, wheat middlings, dehydrated alfalfa meal, soybean oil, corn gluten meal, calcium carbonate, dicalcium phosphate, brewers dried yeast, iodized salt, choline chloride, kaolin, magnesium oxide, L-lysine, DL-methionine, ferrous sulfate, menadione sodium bisulfite complex (source of vitamin K activity), vitamin E acetate, thiamin mononitrate, calcium pantothenate, manganous oxide, niacin, copper sulfate, zinc oxide, vitamin A acetate, pyridoxine hydrochloride, riboflavin, vitamin D₃ supplement, vitamin B₁₂ supplement, folic acid, biotin, calcium iodate, cobalt carbonate.

Macronutrients		
Crude Protein	%	19.1
Fat (ether extract) ^a	%	5.8
Carbohydrate (available) b	%	44.3
Crude Fiber	%	4.6
Neutral Detergent Fiber ^c	%	13.7
Ash	%	6.1
Energy Density ^d	kcal/g (kJ/g)	3.1 (13.0)
Calories from Protein	%	25
Calories from Fat	%	17
Calories from Carbohydrate	%	58
Minerals		
Calcium	%	1.0
Phosphorus	%	0.7
Non-Phytate Phosphorus	%	0.4
Sodium	%	0.3
Potassium	%	0.8
Chloride	%	0.5
Magnesium	%	0.2
Zinc	mg/kg	63
Manganese	mg/kg	93
Copper	mg/kg	23
lodine	mg/kg	3
Iron	mg/kg	240
Selenium	mg/kg	0.16
Amino Acids		
Aspartic Acid	%	1.8
Glutamic Acid	%	2.8
Alanine	%	1.0
Glycine	%	0.8
Threonine	%	0.8
Proline	%	1.4
Serine	%	1.3
Leucine	%	1.7
Isoleucine	%	0.8
Valine	%	0.9
Phenylalanine	%	0.9
Tyrosine	%	0.8
Methionine	%	0.4
Cystine	%	0.3
Lysine	%	1.0
Histidine	%	0.5
Arginine	%	1.2
Tryptophan	%	0.3

Teklad Diets are designed and manufactured for research purposes only.



Standard Product Form: Pellet

- Ctarradra 1	Toddot Form. Fond	
Vitamins		
Vitamin A ^{e, f}	IU/g	30.0
Vitamin D ₃ e, g	IU/g	2.4
Vitamin E	IU/kg	150
Vitamin K ₃ (menadione)	mg/kg	80
Vitamin B ₁ (thiamin)	mg/kg	95
Vitamin B ₂ (riboflavin)	mg/kg	14
Niacin (nicotinic acid)	mg/kg	100
Vitamin B ₆ (pyridoxine)	mg/kg	17
Pantothenic Acid	mg/kg	87
Vitamin B ₁₂ (cyanocobalamin)	mg/kg	0.09
Biotin	mg/kg	0.77
Folate	mg/kg	7
Choline	mg/kg	2200
Fatty Acids		
C16:0 Palmitic	%	0.6
C18:0 Stearic	%	0.2
C18:1ω9 Oleic	%	1.3
C18:2ω6 Linoleic	%	2.6
C18:3ω3 Linolenic	%	0.3
Total Saturated	%	0.8
Total Monounsaturated	%	1.3
Total Polyunsaturated	%	2.9
Other		
Cholesterol	mg/kg	

^a Ether extract is used to measure fat in pelleted diets, while an acid hydrolysis method is required to recover fat in extruded diets. Compared to ether extract, the fat value for acid hydrolysis will be approximately 1% point higher.

For nutrients not listed, insufficient data is available to quantify.

Nutrient data represent the best information available, calculated from published values and direct analytical testing of raw materials and finished product. Nutrient values may vary due to the natural variations in the ingredients, analysis, and effects of processing.

^b Carbohydrate (available) is calculated by subtracting neutral detergent fiber from total carbohydrates.

^c Neutral detergent fiber is an estimate of insoluble fiber, including cellulose, hemicellulose, and lignin. Crude fiber methodology underestimates total fiber.

d Energy density is a calculated estimate of *metabolizable energy* based on the Atwater factors assigning 4 kcal/g to protein, 9 kcal/g to fat, and 4 kcal/g to available carbohydrate.

e Indicates added amount but does not account for contribution from other ingredients.

f 1 IU vitamin A = 0.3 μg retinol

^g 1 IU vitamin D = 25 ng cholecalciferol